

The Difference is Clear.

Water is not enough.

Bench testing shows that water alone **does not** completely remove residue from tube walls, which impedes flow of required nutrients and medications.

THE TUBECLEAR® SYSTEM IS

145.6%

**MORE EFFECTIVE
THAN WATER ALONE!**

Partially Clogged Feeding Tube Treated with Water Flush Only



BEFORE



AFTER

Partially Clogged Feeding Tube Treated with the TubeClear System



BEFORE



AFTER

* Clog contents in the feeding tube were simulated with a feeding formula and fiber supplement and heated to ensure occlusion consistencies between tubes.



Enteral Access Clearing System



When should the TubeClear System be used?

When the feeding tube is sluggish, the nutrition and medication delivery is slowed. Use the TubeClear System to maintain optimal flow in your patients' feeding tubes.



At regular intervals

Bench top results show the TubeClear System is more effective at keeping tube walls clear than water alone.¹ To prevent buildup, use the TubeClear System proactively. This could be after each feeding or once per week.



When tubes become sluggish

Overall, Intensive Care Unit (ICU) patients with feeding tubes receive about 60% of their required nutrition.² When the tube flow is sluggish, use the TubeClear System to restore optimum nutrition and medication delivery.



When tubes become fully blocked

About 25% of all feeding tubes will become clogged.³ Clear the clog with the TubeClear System at bedside to reduce the need to replace the feeding tube, and quickly resume medication and feeding schedules.

1. Actuated Medical, Inc. Internal Test Report Doc. No. 1100791569-000.

2. Guenter, P.; Read, J. "ASPEN Enteral Nutrition by the Numbers: EN Data Across the Healthcare Continuum." American Society for Parenteral and Enteral Nutrition. 2017:20.

3. Dandele, Lauren M. and Lodolee, Amy E. "Efficacy of Agents to Prevent and Treat Enteral Feeding Tube Clogs." The Annals of Pharmacotherapy. 2011; 45:676-680.



Learn more at
TubeClear.com



CONTACT SALES

Actuated Medical, Inc.
+1 (814) 355-0003 x117
info@tubeclear.com
TubeClear.com

Rx ONLY

